




## **APERIO® Hybrid** Thrombectomy Device



Effective hybrid cell design for fast flow restoration  
Outstanding visibility for maximum control and safety  
Broad range of sizes for tailored treatment options

# ORDERING INFORMATION | APERIO® Hybrid

Labelled APERIO® Hybrid Dimensions (mm)	Reference Number	Device Diameter (mm)	Device Length* (mm)	Recommended Vessel Diameter (mm)	Required / Recommended Microcatheters for Delivery (Inch)
3.5 × 28	01-000704	3.5	28	1.5 – 3.0	0.021 NeuroSlider® 21 DLC
4.5 × 30	01-000705	4.5	30	2.0 – 4.0	
4.5 × 40	01-000706	4.5	40	2.0 – 4.0	
4.5 × 50	01-000707	4.5	50	2.0 – 4.0	
6.0 × 40	01-000708	6.0	40	3.5 – 5.5	0.021 – 0.027 NeuroSlider® 21 DLC NeuroSlider® 27 (DLC)
6.0 × 50	01-000709	6.0	50	3.5 – 5.5	

\* Average length within intended vessel diameter

All changes or modifications, may they be technical or other, or changes in the availability of products are expressively reserved.

