



DERIVO® & DERIVO® mini
Embolisation Device



Excellent visibility and perfect adaptability
For vessel diameters from 1.5 mm to 6.0 mm
DERIVO® mini: Up to 3.5 mm device diameter
deliverable through 0.021" ID microcatheters

ORDERING INFORMATION | DERIVO® mini

Labelled DERIVO® mini Dimensions (mm)	Reference Number	Device Diameter (mm)	Device Length (mm)	Recommended Vessel Diameter (mm)	Required / Recommended Microcatheters for Delivery (Inch)
2.5 × 15	01-000428	2.5	15	1.5 – 2.5	0.021 NeuroSlider® 21 NeuroSlider® 21 DLC
2.5 × 20	01-000429	2.5	20	1.5 – 2.5	
2.5 × 20 w/o tip	01-000432	2.5	20	1.5 – 2.5	
2.5 × 25	01-000430	2.5	25	1.5 – 2.5	
2.5 × 25 w/o tip	01-000433	2.5	25	1.5 – 2.5	
3.0 × 15	01-000422	3.0	15	2.0 – 3.0	
3.0 × 20	01-000423	3.0	20	2.0 – 3.0	
3.0 × 20 w/o tip	01-000426	3.0	20	2.0 – 3.0	
3.0 × 25	01-000424	3.0	25	2.0 – 3.0	
3.0 × 25 w/o tip	01-000427	3.0	25	2.0 – 3.0	
3.5 × 15	01-000416	3.5	15	2.5 – 3.5	
3.5 × 20	01-000417	3.5	20	2.5 – 3.5	
3.5 × 20 w/o tip	01-000420	3.5	20	2.5 – 3.5	
3.5 × 25	01-000418	3.5	25	2.5 – 3.5	
3.5 × 25 w/o tip	01-000421	3.5	25	2.5 – 3.5	

ORDERING INFORMATION | DERIVO®

Labelled DERIVO® Dimensions (mm)	Reference Number	Device Diameter (mm)	Device Length (mm)	Recommended Vessel Diameter (mm)	Required / Recommended Microcatheters for Delivery (Inch)
3.5 × 15	01-000408	3.5	15	2.5 – 3.5	0.027 NeuroSlider® 27 (DLC)
3.5 × 20	01-000409	3.5	20	2.5 – 3.5	
3.5 × 25	01-000410	3.5	25	2.5 – 3.5	
3.5 × 30	01-000411	3.5	30	2.5 – 3.5	
3.5 × 40 w/o tip	01-000415	3.5	40	2.5 – 3.5	
4.0 × 15	01-000381	4.0	15	3.0 – 4.0	
4.0 × 20	01-000330	4.0	20	3.0 – 4.0	
4.0 × 25	01-000335	4.0	25	3.0 – 4.0	
4.0 × 30	01-000340	4.0	30	3.0 – 4.0	
4.0 × 40 w/o tip	01-000360	4.0	40	3.0 – 4.0	
4.5 × 15	01-000382	4.5	15	3.5 – 4.5	
4.5 × 20	01-000331	4.5	20	3.5 – 4.5	
4.5 × 25	01-000336	4.5	25	3.5 – 4.5	
4.5 × 30	01-000341	4.5	30	3.5 – 4.5	
4.5 × 40 w/o tip	01-000361	4.5	40	3.5 – 4.5	
5.0 × 15	01-000383	5.0	15	4.0 – 5.0	
5.0 × 20	01-000332	5.0	20	4.0 – 5.0	
5.0 × 25	01-000337	5.0	25	4.0 – 5.0	
5.0 × 30	01-000342	5.0	30	4.0 – 5.0	
5.0 × 40 w/o tip	01-000362	5.0	40	4.0 – 5.0	
5.0 × 50 w/o tip	01-000363	5.0	50	4.0 – 5.0	
5.5 × 15	01-000384	5.5	15	4.5 – 5.5	
5.5 × 20	01-000333	5.5	20	4.5 – 5.5	
5.5 × 25	01-000338	5.5	25	4.5 – 5.5	
5.5 × 30	01-000343	5.5	30	4.5 – 5.5	
5.5 × 40 w/o tip	01-000364	5.5	40	4.5 – 5.5	
5.5 × 50 w/o tip	01-000365	5.5	50	4.5 – 5.5	
6.0 × 15	01-000385	6.0	15	5.0 – 6.0	
6.0 × 20	01-000334	6.0	20	5.0 – 6.0	
6.0 × 25	01-000339	6.0	25	5.0 – 6.0	
6.0 × 30	01-000344	6.0	30	5.0 – 6.0	
6.0 × 40 w/o tip	01-000366	6.0	40	5.0 – 6.0	
6.0 × 50 w/o tip	01-000367	6.0	50	5.0 – 6.0	

Indicated on package as „w/o tip“ as the tip always stays inside the stent for the 40 mm and 50 mm length depending on the vessel diameter.

All changes or modifications, may they be technical or other, or changes in the availability of products are expressly reserved.
Please contact your local Acandis® representative for information on compatible microcatheters.

CE 0297